

Zynex DynaComp

Cold Compression Device

Zynex DynaComp

The DynaComp cold compression device offers targeted relief on-the-go, ensuring patients recover faster and get back to doing what they love.

The DynaComp combines the power of cold therapy and compression in one portable device, making it the perfect solution for quick, effective recovery from work injuries, post-op pain, and more.

Device Features

- ✓ Five Pressure Levels (20-100mmHg)
- ✓ 12 Time Settings (10-120 Minutes)
- ✓ Built-in battery supports 8 hours of use
- ✓ Gel therapy system
- ✓ Small, lightweight, & portable



Treatment Wraps

Wraps are available for:

- Lower Back
- Hip
- Knee
- Shoulder
- Ankle



Two Conservative Treatments

The NexWave electrotherapy device and the DynaComp make an ideal combination for pain management. Combining the two provides the benefits of electrotherapy, cold therapy, and targeted compression at the same time or separately.

