

Zynex Pro OA Knee Brace

Superior Relief for Knee Pain

The Zynex Pro OA Knee Brace (Osteoarthritis) is universally sized— simply choose left or right—to fit virtually any patient. Incredibly easy to apply for both providers and patients. The single upright hinge design is comfortable, lightweight, and functional.

The Zynex Pro OA Knee Brace offloads the medial or lateral compartment for patients with mild to moderate unicompartmental osteoarthritis of the knee. The universal design limits inventory needed to cover the patient size range and allows for the brace to be used for both medial and lateral osteoarthritis.

This OA Knee Brace includes grip padding, which minimizes brace migration or rotation and dramatically increases comfort, wearability, and patient compliance.



Features

- Universal sizing & side specific (left or right)
- Offloads both medial and lateral compartments
- Low profile single upright design
- Precision unloading lock and crank mechanism
- 4-point numbered strapping system

Indications for Use

- Alternative to knee joint replacement
- Degenerative meniscal tears
- Moderate/mild unicompartmental osteoarthritis
- Osteochondral defects
- Articular cartilage defects/repairs



Two Conservative Treatments

Combining knee bracing with electrotherapy offers a supportive approach to managing knee discomfort and maintaining mobility. Knee braces provide stability and support, while electrotherapy may assist with muscle re-education, relaxation of muscle spasms, and increasing local blood circulation. Together, these modalities can help support functional outcomes as part of a broader rehabilitation plan.



Form 30123 Rev 4

Zynex Pro Post-Op Knee Brace

Superior Protection & Stability

The Zynex Pro Post-Op Knee Brace is designed to provide protected, controlled range-of-motion (ROM) for patients recovering from knee surgery or those who have knee injuries or instabilities.

Streamlined Design

- Telescoping struts for length adjustment from 17" - 29" (43.18 cm - 76.33 cm)
- Simple rotate and release flexion/extension settings
- Unique safety lock to increase patient compliance
- 100% tool-free adjustments
- Quickly locks at 0°, 15°, 30°, 45°

Sizing Options

Standard - Up to 29" (76.33 cm)
X-Large - Up to 35" (88.9 cm)



Indications for Use

- Post-op rehabilitation
- Soft tissue knee injuries
- Knee sprains/strains
- Tibial plateau fractures
- Osteochondral repairs
- Meniscal repairs
- Patella tendon repairs
- Condylar fractures
- ACL, PCL, MCL, LCL repairs/injuries
- HTOs (High Tibial Osteotomy)
- Stable or internally fixed fractures of the knee

Two Conservative Treatments

Combining knee bracing with electrotherapy offers a supportive approach to managing knee discomfort and maintaining mobility. Knee braces provide stability and support, while electrotherapy may assist with muscle re-education, relaxation of muscle spasms, and increasing local blood circulation. Together, these modalities can help support functional outcomes as part of a broader rehabilitation plan.

