

# Zynex Pro Post-Op Knee Brace

## Superior Protection & Stability

The Zynex Pro Post-Op Knee Brace is designed to provide protected, controlled range-of-motion (ROM) for patients recovering from knee surgery or those who have knee injuries or instabilities.

## Streamlined Design

- Telescoping struts for length adjustment from 17" - 29" (43.18 cm - 76.33 cm)
- Simple rotate and release flexion/extension settings
- Unique safety lock to increase patient compliance
- 100% tool-free adjustments
- Quickly locks at 0°, 15°, 30°, 45°

## Sizing Options

Standard - Up to 29" (76.33 cm)  
X-Large - Up to 35" (88.9 cm)



## Indications for Use

- Post-op rehabilitation
- Soft tissue knee injuries
- Knee sprains/strains
- Tibial plateau fractures
- Osteochondral repairs
- Meniscal repairs
- Patella tendon repairs
- Condylar fractures
- ACL, PCL, MCL, LCL repairs/injuries
- HTOs (High Tibial Osteotomy)
- Stable or internally fixed fractures of the knee

## Two Conservative Treatments

Combining knee bracing with electrotherapy offers a supportive approach to managing knee discomfort and maintaining mobility. Knee braces provide stability and support, while electrotherapy may assist with muscle re-education, relaxation of muscle spasms, and increasing local blood circulation. Together, these modalities can help support functional outcomes as part of a broader rehabilitation plan.

