

Zynex TLSO

Thoracic-Lumbar-Sacral Orthosis

Zynex TLSO

The Zynex TLSO is a dual-purpose orthosis that supports, stabilizes, and corrects the mid and lower spine (T9-S1).

The Zynex TLSO's dual pulley system telescopes the thoracic core independently from the lumbar brace. This revolutionary mechanism secures a neutral spine alignment, and an increased low load assisted thoracic stretch.

Additionally, the Zynex TLSO can be used to provide unmatched rotational stability and motion restriction during the post-operative period.

Key Features

- Effectively assist postural correction and reeducation
- Provides rotational stability and motion restriction
- Low profile design
- No grip closure for patients lacking hand strength
- Adjustable chest strap for a custom fit
- Patient-friendly design



Indications for Use



Spinal Stenosis



Herniated Disc



Diagnosis requiring immobilization of the thoracic & lumbar spine



Degenerative Spinal Fractures



Kyphosis

Size	Waist Meas.	Size	Waist Meas.
0	22" - 32"	4	42" - 46"
1	31" - 37"	5	45" - 49"
2	36" - 40"	6	48" - 52"
3	39" - 43"	7	51" - 55"

Sizing Guide

- Universal (XS - 3XL) / (22" - 55")
- 10" Extender panel available for patients needing an excess of 55"

